


SARS-CoV-2 HUMAN IgG QUANTITATIVE (4-PLEX)



The SARS-CoV-2 Human IgG Quantitative (4-Plex) array is a multiplexed serological assay to be used in research and larger surveillance studies.

Q-Plex™ Human IgG Quantitative (4-Plex) Array

Key Attributes:

Identifies individuals with an adaptive immune response to SARS-CoV-2, indicating prior infection

97% sensitivity, in comparison to the gold standard of PCR

Assay recognizes IgG antibodies to subunits S1 and S2 of the viral Spike protein

Reactivity to both S1 and S2 reduces false negative results, increasing assay sensitivity

Results in 2 hours

ELISA based assay

96-well format compatible with automation

The Q-Plex SARS-CoV-2 Human IgG (4-plex) Assay is available for immediate use

Test Applications

Vaccine development- quantify antibody response to spike protein subunits

Determine antibody levels post COVID-19 recovery



Determine COVID-19 prevalence in a defined population

Quantify individual antibody response to an infection or a vaccination

General Information

The Q-Plex™ SARS-CoV-2 Human IgG (4-plex) array simultaneously recognizes human IgG antibodies to S1 and S2 of the spike protein. The spike protein recognizes a cell surface receptor to which the virus binds and inserts its “coding” information into the cell. Inclusion of the S2 subunit as a target, in this assay, decreases the false positive test rate, providing additional test accuracy. “High quality tests are essential as we work to better understand the SARS-CoV-2 pandemic. In addition, improved quality and accuracy of the assay is vital to correctly evaluate the spread of this pandemic, which varies across communities,” noted Adam Brown, CEO of Quansys Biosciences.

Performance Characteristics

SARS-CoV-2 IgG Assay vs. Molecular COVID-19 Test		
N=254	Confirmed Positive	Confirmed Negative
IgG Test Positive	33	0
IgG Test Negative	1	220

Estimated Sensitivity (PPA)	97%
Estimated Specificity (NPA)	100%
PPV	100%
NPV	100%

Note: A sample is considered positive for reactive antibodies to SARS-CoV-2 when both the S1 and S2 concentrations are above their respective cutoffs.

Human SARS-CoV-2 S1 IgG Quantitative Assay

Assay Range:	1000 - 1.3 U/mL
Lower Limit of Detection:	1.3 U/mL
Negative/Positive Cutoff:	7 U/ml
Precision:	Intra-assay 9%, Inter-assay 13%
Linearity:	93-113%, average 102%

Human SARS-CoV-2 S2 IgG Quantitative Assay

Assay Range:	1000 - 1.3 U/mL
Lower Limit of Detection:	1.3 U/mL
Negative/Positive Cutoff:	24 U/ml
Precision:	Intra-assay 8%, Inter-assay 10%
Linearity:	92-107%, average 100%

Precision

Intra-assay (n=18)			Inter-assay (n=6)		
	Concentration (U/mL)	%CV		Concentration (U/mL)	%CV
Sample 1	6.5	12%	Sample 1	106.75	10.2%
Sample 2	118.2	7%	Sample 2	6.36	15.4%
Sample 3	1.1	6%	Sample 3	5.49	13.0%
Average Intra-assay CV:		9%	Average Inter-assay CV:		13%

Intra-assay (n=18)			Inter-assay (n=6)		
	Concentration (U/mL)	%CV		Concentration (U/mL)	%CV
Sample 1	6.5	12%	Sample 1	106.75	10.2%
Sample 2	118.2	7%	Sample 2	6.36	15.4%
Sample 3	1.1	6%	Sample 3	5.49	13.0%
Average Intra-assay CV:		9%	Average Inter-assay CV:		13%

Dilutional Linearity

Dilution Factor	% Recovery		
	1:100	1:200	1:400
Sample 1	97%	100%	97%
Sample 2	100%	111%	113%
Sample 3	108%	101%	93%
Average % Linearity:			102%
% Linearity Range:			93 - 113%

Dilution Factor	% Recovery		
	1:100	1:200	1:400
Sample 1	104%	93%	95%
Sample 2	97%	106%	106%
Sample 3	107%	102%	92%
Average % Linearity:			100%
% Linearity Range:			92 - 107%

SARS-CoV-2 Kit

Catalog #	Product Name	
711649HU	SARS-CoV-2 Human IgG Quantitative (4-Plex)	Call for pricing

SARS-CoV-2 Human IgG (4-Plex) details	
Species:	Human
Assay Type:	Indirect
Detection Method:	Chemiluminescent
Sample Type:	Human Serum, Plasma, Whole Blood
Assay Length:	2 hours, RT
Volume Required Per Well:	Min 2µl
Total Wash Steps:	3
Multiplex Format:	96-well solid plate
Catalog No:	711649HU




Q-View™ Imager generated image demonstrating typical sample results and reproducibility (samples shown were run in duplicates in each row). Samples can be run in duplicate or singlet allowing 40 to 80 unique samples to be run simultaneously on the same plate.



Why IgG

With the goal of producing a highly specific and precise assay, Quansys' R&D intentionally developed the SARS-CoV-2 antibody test to detect IgGs from non-acute disease phase individuals. IgM testing serves a purpose in understanding the initial immune response days after and during infection. However, known issues with IgM testing specificity coupled with the natural decline of IgM levels through disease course diminishes their utility for periods beyond the first few weeks of infection. Focusing on IgGs, simplifies interpretation, and more easily permits the flexibility of the Q-Plex platform to add additional viral antigens from the parent strain or mutant strains, while maintaining high levels of specificity. Moreover, the Q-Plex platform is built for quantitation and long term/periodic monitoring of IgG levels will be possible when calibration standards are established.



flexibility of the Q-Plex platform to add additional viral antigens from the parent strain or mutant strains, while maintaining high levels of specificity. Moreover, the Q-Plex platform is built for quantitation and long term/periodic monitoring of IgG levels will be possible when calibration standards are established.



Quansys Biosciences
 365 North 600 West
 Logan, UT 84321
www.quansysbio.com
 Fax: 435-750-6869
 Phone: 435-752-0531
 Toll Free: 888-QUANSYS (782-6797)



Quansys is an ISO 9001:2015 and ISO 13485:2016 registered company and complies with cGMP. Products are designed, developed, and manufactured according to the procedures outlined in our Quality Management System.

©2020 Quansys Biosciences Inc. All rights reserved. Q-View and Q-Plex are Trademarks of Quansys Biosciences Inc. ISO is a trademark of the International Standards Organization. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

