

# Q-PLEX™ MOUSE CYTOKINE SCREEN HS STRIPWELLS (15-PLEX)



The Q-Plex™ Mouse Cytokine HS Stripwells (15-Plex) is a fully quantitative ELISA-based chemiluminescent assay allowing the concurrent measurement of sixteen biomarkers or analytes.

**QUANSYS**   
B I O S C I E N C E S

Catalog #	Product	GMCSF	IFN $\gamma$	IL-1 $\alpha$	IL-1 $\beta$	IL-2
-----------	---------	-------	--------------	---------------	--------------	------

112149MS-15P	Q-Plex™ Mouse Cytokine Screen HS Stripwells (15-Plex)	IL-3	IL-4	IL-5	IL-6	IL-10
		IL-12p70	IL-17A	MIP-1 $\alpha$	RANTES	TNF $\alpha$

# Background

The Q-Plex™ Mouse Cytokine HS Stripwells (15-Plex) is a fully quantitative ELISA-based chemiluminescent assay allowing the concurrent measurement of sixteen biomarkers or analytes: GMCSF, IFN $\gamma$ , IL-1 $\alpha$ , IL-1 $\beta$ , IL-2, IL-3, IL-4, IL-5, IL-6, IL-10, IL-12p70, IL-17A, MIP-1 $\alpha$ , RANTES, and TNF $\alpha$ .

Each kit contains a 96-well plate, featuring the relevant biomarker panel, and all the reagents required to perform testing. Q-Plex™ plates are built by absorbing sixteen distinct capture antibodies in a defined array to the bottom of each well. Our high quality reagents help ensure the accuracy of your results.

# Example Performance and Validation Data

Analyte	Assay Type	Calibrator Range	Upper Limit of Quantification (ULOQ)	Lower Limit of Quantification (LLOQ)	Limit of Detection	Precision* (Inter-assay)	Precision* (Intra-assay)	Average* Linearity
GMCSF	Sandwich	560 - 0.77 (pg/mL)	392 (pg/mL)	1 (pg/mL)	0.38 (pg/mL)	4.9%	7.5%	94%
IFN $\gamma$	Sandwich	516 - 0.71 (pg/mL)	361 (pg/mL)	0.94 (pg/mL)	0.35 (pg/mL)	3.3%	5.5%	110%
IL-1 $\alpha$	Sandwich	610 - 0.84 (pg/mL)	427 (pg/mL)	1.1 (pg/mL)	0.42 (pg/mL)	3.3%	7.4%	108%
IL-1 $\beta$	Sandwich	2500 - 3.43 (pg/mL)	1750 (pg/mL)	4.5 (pg/mL)	1.71 (pg/mL)	3.1%	9.3%	81%
IL-2	Sandwich	360 - 0.49 (pg/mL)	252 (pg/mL)	0.6 (pg/mL)	0.25 (pg/mL)	5.1%	6.6%	88%
IL-3	Sandwich	60 - 0.08 (pg/mL)	42 (pg/mL)	0.1 (pg/mL)	0.04 (pg/mL)	4.7%	7.0%	106%
IL-4	Sandwich	400 - 0.55 (pg/mL)	280 (pg/mL)	0.7 (pg/mL)	0.27 (pg/mL)	6.2%	6.6%	98%
IL-5	Sandwich	546 - 0.75 (pg/mL)	382 (pg/mL)	0.96 (pg/mL)	0.37 (pg/mL)	6.6%	5.2%	103%
IL-6	Sandwich	504 - 0.69 (pg/mL)	353 (pg/mL)	0.91 (pg/mL)	0.35 (pg/mL)	7.3%	6.1%	106%
IL-10	Sandwich	667 - 0.91 (pg/mL)	467 (pg/mL)	1.33 (pg/mL)	0.46 (pg/mL)	7.5%	3.0%	108%
IL-12p70	Sandwich	572 - 0.78 (pg/mL)	400 (pg/mL)	1.01 (pg/mL)	0.39 (pg/mL)	5.5%	10.8%	107%
IL-17A	Sandwich	760 - 1.04 (pg/mL)	532 (pg/mL)	1.4 (pg/mL)	0.52 (pg/mL)	6.5%	6.1%	92%
MIP-1 $\alpha$	Sandwich	1000 - 1.37 (pg/mL)	700 (pg/mL)	1.8 (pg/mL)	0.69 (pg/mL)	10.9%	6.7%	86%
RANTES	Sandwich	1653 - 2.27 (pg/mL)	1157 (pg/mL)	2.92 (pg/mL)	1.13 (pg/mL)	5.0%	5.5%	108%
TNF $\alpha$	Sandwich	429 - 0.59 (pg/mL)	300 (pg/mL)	0.69 (pg/mL)	0.29 (pg/mL)	6.9%	5.2%	104%

\* Actual values may vary from kit to kit. Please see the In-Kit Certificate of Analysis included in your kit for specific values.



### Quansys Biosciences

365 North 600 West  
 Logan, UT 84321  
[www.quansysbio.com](http://www.quansysbio.com)

Fax: 435-750-6869  
 Phone: 435-752-0531  
 Toll Free: 888-QUANSYS (782-6797)



Quansys is an ISO 9001:2015 and ISO 13485:2016 registered company and complies with GMP. Products are designed, developed, and manufactured according to the procedures outlined in our Quality Management System.

©2024 Quansys Biosciences Inc. All rights reserved. Q-View and Q-Plex are Trademarks of Quansys Biosciences Inc. ISO is a trademark of the International Standards Organization. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.