

## CERTIFICATE OF ANALYSIS

Product Name: Q-Plex™ Human  
Angiogenesis (9-Plex)


Order No. 21575

**REF** Part No. 150233HU

**LOT** Lot No. HAGK200518 17735 001

**Exp Date:** 2023-MAR-21

**Software Product Code:** HA1M200518-AGQ



### Reconstitute Calibrator with 650 µL of Sample Diluent

Component	<b>REF</b> Part No.	<b>LOT</b> Lot No.	<b>Exp. Date</b>
Plate	150252HU	HAGP220222	2024-FEB-22
Sample Diluent	101104HU	HQ6S210519	2023-MAY-19
Calibrator	102655HU	HA1M200518	2024-APR-12
Competitor	N/A	N/A	N/A
Detection	150256HU_01	HAGDQ6220302	2024-MAR-02
Streptavidin HRP	101173GR	QNSN210921	2023-MAR-21
Substrate A	101159GR	QNSA220310	2024-MAR-10
Substrate B+	101120GR	QNSB+220310	2024-MAR-10
Wash Buffer	101158GR	QNSW220120	2024-JAN-20

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	ANG-2	pg/mL	8600	8600	11.80	5.90
2	FGF	pg/mL	5200	5200	7.13	3.57
3	HGF	pg/mL	13700	13700	18.79	9.40
4	IL-8	pg/mL	1800	1800	2.47	1.23
5	PDGF	pg/mL	4300	4300	5.90	2.95
6	TIMP-1	pg/mL	14193	14200	19.48	9.74
7	TIMP-2	pg/mL	18200	18200	24.97	12.48
8	TNFa	pg/mL	5400	5400	7.41	3.70
9	VEGF	pg/mL	2000	2000	2.74	1.37
10	REF	N/A	N/A	N/A	N/A	N/A

Quality Release:

*Stacey Pridmore*

101896GR 2.0

for Research Use Only