

CERTIFICATE OF ANALYSIS

Product Name: Q-Plex™ Human
Chemokine (9-Plex)


Order No. 21624

REF Part No. 120233HU

LOT Lot No. HCKK200518 17772 001

 Exp Date: 2023-MAR-01

Software Product Code: HA1M200518-CKQ



Reconstitute Calibrator with 650 µL of Sample Diluent

Component	REF Part No.	LOT Lot No.	 Exp. Date
Plate	140252HU	HCKP220202	2024-FEB-02
Sample Diluent	101104HU	HQ6S210301	2023-MAR-01
Calibrator	102655HU	HA1M200518	2024-APR-21
Competitor	N/A	N/A	N/A
Detection	120256HU_01	HCKDQ6210804	2023-AUG-04
Streptavidin HRP	101173GR	QNSN210921	2023-MAR-21
Substrate A	101159GR	QNSA220309	2024-MAR-09
Substrate B+	101120GR	QNSB+220309	2024-MAR-09
Wash Buffer	101158GR	QNSW220120	2024-JAN-20

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	Eotaxin	pg/mL	1600	1600	2.19	1.10
2	GROa	pg/mL	2000	2000	2.74	1.37
3	I-309	pg/mL	2400	2400	3.29	1.65
4	IL-8	pg/mL	1800	1800	2.47	1.23
5	IP-10	pg/mL	3000	1000	4.12	2.06
6	MCP-1	pg/mL	1900	633	2.61	1.30
7	MCP-2	pg/mL	510	510	6.30	3.15
8	RANTES	pg/mL	6200	6200	8.50	4.25
9	TARC	pg/mL	1400	1400	1.92	0.96
10	REF	N/A	N/A	N/A	N/A	N/A

Quality Release:

Stacey Priober
for Research Use Only

101896GR 2.0