

## CERTIFICATE OF ANALYSIS

Product Name: Q-Plex™ Human TGFb HDR


Order No. 21575

**REF** Part No. 447249HU


**LOT** Lot No. HTGBK181029 17735 002

 Exp Date: 2023-MAR-21

**Software Product Code: HTGBM181029TG**



**Reconstitute Calibrator with 550 µL of Sample Diluent**

Component	<b>REF</b> Part No.	<b>LOT</b> Lot No.	 Exp. Date
Plate	447252HU	HTGBP210512	2023-MAY-12
Sample Diluent	101104HU	HQ6S210519	2023-MAY-19
Calibrator	447255HU	HTGBM181029	2024-APR-12
Competitor	N/A	N/A	N/A
Detection	447256HU	HTGBD210525	2023-MAY-25
Streptavidin HRP	101173GR	QNSN210921	2023-MAR-21
Substrate A	101159GR	QNSA220310	2024-MAR-10
Substrate B+	101120GR	QNSB+220310	2024-MAR-10
Wash Buffer	101158GR	QNSW220120	2024-JAN-20

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	TGFb	pg/mL	31100	31100	1.99	1.00
2	TGFb	pg/mL	31100	31100	1.99	1.00
3	TGFb	pg/mL	31100	31100	1.99	1.00
4	TGFb	pg/mL	31100	31100	1.99	1.00
5	TGFb	pg/mL	31100	31100	1.99	1.00
6	TGFb	pg/mL	31100	31100	1.99	1.00
7	REF	N/A	N/A	N/A	N/A	N/A

Quality Release: \_\_\_\_\_

*Stacey Prodan*

101896GR 2.0

for Research Use Only