

CERTIFICATE OF ANALYSIS

Product Name: Q-Plex™ Human Custom (2-Plex)


Order No. 21623

REF Part No. 107749GR

LOT Lot No. HCUK220526 17908 001

Exp Date: 2023-MAR-01

Software Product Code: **HCUM220527-AP**



Reconstitute Calibrator 1 with 750 µL of Complete Sample Diluent

Reconstitute Calibrator 2 with 300 µL of Calibrator 1

Component	REF Part No.	LOT Lot No.	Exp. Date
Plate	101852GR	HCUP220427	2024-APR-27
Sample Diluent	101104HU	HQ6S210301	2023-MAR-01
Calibrator 2	101855GR	HCUM220526	2024-MAY-31
Calibrator 1	112455HU	HCHM200803	2024-MAY-24
Detection	101856GR	HCUD220527	2024-MAY-27
Streptavidin HRP	101173GR	QNSN220111	2023-JUL-11
Substrate A	101159GR	QNSA220309	2024-MAR-09
Substrate B+	101120GR	QNSB+220309	2024-MAR-09
Wash Buffer	101158GR	QNSW211201	2023-DEC-01

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	CD-163	ng/mL	805	805	3.31276	1.65638
2	IL-6 HS	pg/mL	240	240	80	40
3	REF	N/A	N/A	N/A	N/A	N/A

Quality Release:

Stephen Adams
for Research Use Only

101896GR 2.0