

Q-Plex™ ARRAY

Q-Plex™ SARS-CoV-2 Human IgG (4-Plex) Kit

Software Product Code: HCVM200410


S1 Cutoff Value: 1.3

S2 Cutoff Value: 0.4

REF	691649HU
LOT	HCVK200410 17830 001
	2022-DEC-14

Product Information				
Component	Part Number	Lot Number	Exp. Date	Recon. Volume
Plate	691652HU	HCVP220106	2024-JAN-06	N/A
Sample Diluent	691657HU	HCVS211214	2022-DEC-14	N/A
High Control	691631HU	HCVM200410-H	2023-NOV-09	250 µL
Low Control	691632HU	HCVM200410-L	2023-NOV-09	250 µL
Negative Control	691643HU	HCVM200410-N	2023-NOV-09	250 µL
Detection	691643HU	HCVND211210	2024-DEC-10	N/A
Streptavidin HRP	691656HU	QNSN210721	2023-JAN-21	N/A
Substrate A	691173GR	QNSA211206	2023-DEC-06	N/A
Substrate B+	691159GR	QNSB+211206	2023-DEC-06	N/A
Wash Buffer	691120GR	QNSW211201	2023-DEC-01	N/A

Assay Performance			
Analyte	High Control Range (U/mL)	Low Control Range (U/mL)	Negative Control Range (U/mL)
S1	13629	1035	121
S2	16710	6091	964
Sheep Fc	100	100	100
Hu IgG	59747	60622	61627



Quality Control:

Sophia Adams

Date:

17 MAY 22

Quality Assurance:

[Signature]

Date:

17 MAY 22