


CERTIFICATE OF ANALYSIS

Product Name: Q-Plex™ Human
Enviromental Enteric
Dysfunction (11-Plex)

Order No. 19531
Part No. 569849HU
Lot No. HMDK210721 16826 001
Exp Date: 2023-JUN-02

REF
LOT

Software Product Code: HMDM210721-PTH



Reconstitute Calibrator with 300 µL of Sample Diluent

Component	REF Part No.	LOT Lot No.	Exp. Date
Plate	569852HU	HMDP220524	2024-MAY-24
Sample Diluent	101104HU	HQ6S220107	2024-JAN-07
Calibrator	569855HU	HMDM210721	2024-JUN-20
Competitor	569890HU	HMDC210708	2024-JUN-20
Detection	569856HU	HMDD220602	2023-JUN-02
Streptavidin HRP	101173GR	QNSN220321	2023-SEP-21
Substrate A	101159GR	QNSA220202	2024-FEB-02
Substrate B2	101128GR	QNSB2220202	2024-FEB-02
Wash Buffer	101158GR	QNSW211201	2023-DEC-01

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	Ferritin	µg/L	111	111	0.15	0.076
2	Tg	µg/L	12.4	12.4	0.017	0.0085
3	sTfR	mg/L	75	75	0.10	0.0514
4	HRP2	µg/L	0.96	0.96	0.0013	0.00066
5	IGF-1	ng/mL	202	202	0.28	0.14
6	FGF-21	pg/mL	688	688	0.94	0.47
7	IFABP	pg/mL	2145	2145	2.94	1.471
8	REF	N/A	N/A	N/A	N/A	N/A
9	CRP	mg/L	5.37	5.37	0.022	0.011
10	CD-14	ng/mL	2170	2170	8.93	4.47
11	AGP	g/L	0.30	0.30	0.0012	0.00062
12	RBP4	µmol/L	1.53	1.53	0.0063	0.0031

Quality Release: *Hacey Priston*
for Research Use Only

101896GR 2.0