

CERTIFICATE OF ANALYSIS

Product Name: Q-Plex™ Mouse
Custom (4-Plex)


Order No. 21669A

REF Part No. 107749GR

LOT Lot No. MCUK191126 17945 001

Exp Date: 2023-NOV-20

Software Product Code: MA1M191126-CM4



Reconstitute Calibrator with 550 µL of Sample Diluent

Component	REF Part No.	LOT Lot No.	Exp. Date
Plate	101852GR	MCUP220609	2024-JUN-09
Sample Diluent	110957MS	MCSS220607	2024-JUN-07
Calibrator	102655MS	MA1M191126	2024-JUN-16
Competitor	N/A	N/A	N/A
Detection	101856GR	MCUD220610	2024-JUN-10
Streptavidin HRP	101173GR	QNSN220520	2023-NOV-20
Substrate A	101159GR	QNSA211217	2023-DEC-17
Substrate B+	101120GR	QNSB+211217	2023-DEC-17
Wash Buffer	101158GR	QNSW220422	2024-APR-22

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	IL-1b	pg/mL	12950	N/A	N/A	12.41
2	IL-6	pg/mL	4800	N/A	N/A	2.90
3	MCP-1	pg/mL	3150	N/A	N/A	3.40
4	TNFa	pg/mL	2980	N/A	N/A	3.40
5	REF	N/A	N/A	N/A	N/A	N/A

Quality Release:

Sophia Adams
for Research Use Only

101896GR 2.0