

CERTIFICATE OF ANALYSIS

Product Name: Q-Plex™ Human Cytokine Panel
2 (6-Plex)


Order No. 21856

REF Part No. 115333HU

LOT Lot No. HCP3K200518 18028 001

Exp Date: 2023-MAR-01

Software Product Code: HA1M200518-P26P



Reconstitute Calibrator with 650 µL of Sample Diluent

Component	REF Part No.	LOT Lot No.	Exp. Date
Plate	115352HU	HCP3P220217	2024-FEB-17
Sample Diluent	101104HU	HQ6S210301	2023-MAR-01
Calibrator	102655HU	HA1M200518	2024-JUL-01
Competitor	N/A	N/A	N/A
Detection	110956HU_01	HCSDQ6220407	2024-APR-07
Streptavidin HRP	101173GR	QNSN220520	2023-NOV-20
Substrate A	101159GR	QNSA211217	2023-DEC-17
Substrate B+	101120GR	QNSB+211217	2023-DEC-17
Wash Buffer	101158GR	QNSW211201	2023-DEC-01

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	IL-6	pg/mL	3100	3100	4.25	2.13
2	IL-8	pg/mL	1800	1800	2.47	1.23
3	TNFA	pg/mL	5400	5400	7.41	3.70
4	IFNγ	pg/mL	8000	8000	10.97	5.49
5	IL-4	pg/mL	1600	1600	2.19	1.10
6	IL-10	pg/mL	6800	6800	9.33	4.66
7	REF	N/A	N/A	N/A	N/A	N/A

Quality Release:

Stacey Prater

101896GR 2.0

for Research Use Only