


CERTIFICATE OF ANALYSIS

Product Name: Q-Plex™ Human
Custom (8-Plex)
Order No. 21664
Part No. 107749GR
Lot No. HCUK200131 18121 001
Exp Date: 2023-JUN-13
Software
Product Code: FHCM190715-IL



Reconstitute Calibrator 1 with 550 µL of Complete Sample Diluent

Reconstitute Calibrator 2 with 350 µL of Calibrator 1

Component	REF Part No.	LOT Lot No.	Exp. Date
Plate	101852GR	HCUP220613	2023-JUN-13
Sample Diluent	101122HU	HQ16S220628	2023-JUN-28
Calibrator 1	332155HU	FHCM190715	2023-JUL-20
Calibrator 2	101855GR	HCUM200131	2023-JUL-20
Competitor	101890GR	HCUC220616	2023-JUL-20
Detection	101856GR	HCUD220616	2023-JUN-16
Streptavidin HRP	101173GR	QNSN220711	2024-JAN-11
Substrate A	101159GR	QNSA220620	2023-DEC-20
Substrate B+	101120GR	QNSB+220620	2023-DEC-20
Wash Buffer	101158GR	QNSW220615	2023-DEC-15

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	*Cortisol	ng/mL	714	N/A	N/A	N/A
2	C-Peptide	ng/mL	36	N/A	N/A	N/A
3	*CRP	ng/mL	569	N/A	N/A	N/A
4	*E1G	ng/mL	522	N/A	N/A	N/A
5	FSH	mIU/mL	193	N/A	N/A	N/A
6	HCGb	ng/mL	40	N/A	N/A	N/A
7	IL-6	pg/mL	1396	N/A	N/A	N/A
8	IL-10	pg/mL	607	N/A	N/A	N/A
9	REF	N/A	N/A	N/A	N/A	N/A

*Competitive Assay

Quality Release:

101896GR 2.0

for Research Use Only