


CERTIFICATE OF ANALYSIS

Product Name: Q-Plex™ Human Malaria (5-Plex)
 Order No. 21629-2
 Part No. 565949HU
 Lot No. HMRK211108 17931 001
 Exp Date: 2023-MAY-15
 Software Product Code: HMRM211108



Reconstitute Calibrator with 375 µL of Calibrator Diluent

Component	REF	Part No.	LOT	Lot No.	Exp. Date
Plate		565952HU		HMRP220606	2023-DEC-06
Caibrator Diluent		565928HU		HMRU211115	2023-MAY-15
Assay Diluent		565983HU		HMRV211116	2023-MAY-16
Calibrator		565955HU		HMRM211108	2023-JUL-29
Competitor		565990HU		HMRC191014	2023-JUL-29
Detection		565956HU		HMRD211119	2023-MAY-19
Streptavidin HRP		101173GR		QNSN220111	2023-JUL-11
Moss Substrate A+		105959GR		MOSA+220714	2024-JUL-14
Moss Substrate B		105920GR		MOSB220714	2024-JUL-14
Wash Buffer		101158GR		QNSW211201	2023-DEC-01

Spot #	Analyte	Units	Cal Conc	ULOQ	LLQ	LLD
1	HRP2	pg/mL	2800	2800	0.68	0.34
2	LDH-Pan	pg/mL	67000	67000	16.36	8.18
3	LDH-Pv	pg/mL	19200	19200	4.69	2.34
4	LDH-Pf	pg/mL	20800	20800	5.08	2.54
5	*CRP	ng/mL	38000	38000	37.11	18.55
6	REF	N/A	N/A	N/A	N/A	N/A

***Competitive Assay**

Quality Release:

[Signature]

101896GR 2.0

for Research Use Only