


## CERTIFICATE OF ANALYSIS

Product Name: Q-Plex™ Mouse Cytokine  
Panel 2 (4-Plex)


Order No. 21933

**REF** Part No. 115549MS


**LOT** Lot No. MCP4K191126 18111 002

 Exp Date: 2024-JAN-11

Software Product Code: **MA1M191126-P2A**



**Reconstitute Calibrator with 550 µL of Sample Diluent**

Component	<b>REF</b> Part No.	<b>LOT</b> Lot No.	 Exp. Date
Plate	115552MS	MCP4P220322	2024-MAR-22
Sample Diluent	110957MS	MCSS220602	2024-JUN-02
Calibrator	102655MS	MA1M191126	2024-JUL-19
Competitor	N/A	N/A	N/A
Detection	101856GR	MCUD220706	2024-JUL-06
Streptavidin HRP	101173GR	QNSN220711	2024-JAN-11
Substrate A	101159GR	QNSA220620	2024-JUN-20
Substrate B+	101120GR	QNSB+220620	2024-JUN-20
Wash Buffer	101158GR	QNSW220615	2024-JUN-15

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	IFN $\gamma$	pg/mL	7130	N/A	N/A	5.40
2	IL-1b	pg/mL	12950	N/A	N/A	12.41
3	IL-6	pg/mL	4800	N/A	N/A	2.90
4	TNF $\alpha$	pg/mL	2980	N/A	N/A	3.40
5	REF	N/A	N/A	N/A	N/A	N/A

Quality Release:

*Alicia Proctor*

101896GR 2.0

for Research Use Only