

CERTIFICATE OF ANALYSIS

Product Name: Q-Plex™ Porcine
Cytokine - High
Sensitivity (4-Plex)

Order No. 21914

REF

Part No. 119149PC

LOT

Lot No. PCSK170306 18095 001


Exp Date:

2023-JUN-15

Software Product

Code:

PCSM170306CY



Reconstitute Calibrator with 500 µL of Sample Diluent

Component	REF Part No.	LOT Lot No.	Exp. Date
Plate	119152PC	PCSP210726	2023-JUL-26
Sample Diluent	119157PC	PCSS220321	2024-MAR-21
Calibrator	119155PC	PCSM170306	2024-JUL-12
Competitor	N/A	N/A	N/A
Detection	119156PC	PCSD210615	2023-JUN-15
Streptavidin HRP	101173GR	QNSN220520	2023-NOV-20
Substrate A	101159GR	QNSA220202	2024-FEB-02
Substrate B+	101120GR	QNSB+220202	2024-FEB-02
Wash Buffer	101158GR	QNSW220615	2024-JUN-15

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	IL-1b	pg/mL	1600	N/A	N/A	2.10
2	IL-6	pg/mL	2000	N/A	N/A	2.50
3	IL-8	pg/mL	2300	N/A	N/A	2.90
4	TNFa	pg/mL	1350	N/A	N/A	1.90
5	REF	N/A	N/A	N/A	N/A	N/A

Quality Release:

[Signature]
for Research Use Only

101896GR 2.0