


Q-Plex™ ARRAY

CERTIFICATE OF ANALYSIS

Product Name: Q-Plex™ Human Custom
 Singleplex
 Order No. 22012
 REF Part No. 107749GR
 LOT Lot No. HCUK220526 18284 001
 Exp Date: 2023-DEC-01
 Software
 Product Code: **HCUM220526-CD163**



Reconstitute Calibrator with 300 µL of Complete Sample Diluent

Component	Part No.	Lot No.	Exp. Date
Plate	101852GR	HCUP220815	2024-AUG-15
Sample Diluent	101104HU	HQ6S220107	2024-JAN-07
Calibrator	101855GR	HCUM220526	2024-AUG-19
Competitor	101890GR	HCUC220527	2024-AUG-19
Detection	N/A	N/A	N/A
Streptavidin HRP	101173GR	QNSN220715	2024-JAN-15
Substrate A	101159GR	QNSA220713	2024-JUL-13
Substrate B+	101120GR	QNSB+220713	2024-JUL-13
Wash Buffer	101158GR	QNSW211201	2023-DEC-01

Spot #	Analyte	Units	Cal	ULOQ	LLOQ	LLD
1	CD-163	ng/mL	805	805	3.31	1.66
2	CD-163	ng/mL	805	805	3.31	1.66
3	REF	N/A	N/A	N/A	N/A	N/A

Quality Release:  101896GR 2.0
 for Research Use Only