


CERTIFICATE OF ANALYSIS

Product Name: Q-Plex™ Mouse Cytokine
 Panel 2 (4-Plex)
 Order No. 21952
REF Part No. 115549MS
LOT Lot No. MCP4K191126 18135 001
 Exp Date: 2024-JAN-11
 Software
 Product Code: **MA1M191126-P2A**



Reconstitute Calibrator with 550 µL of Sample Diluent

Component	Part No.	Lot No.	Exp. Date
Plate	115552MS	MCP4P220801	2024-AUG-01
Sample Diluent	110957MS	MCSS220602	2024-JUN-02
Calibrator	102655MS	MA1M191126	2024-AUG-05
Competitor	N/A	N/A	N/A
Detection	101856GR	MCUD220706	2024-JUL-06
Streptavidin HRP	101173GR	QNSN220711	2024-JAN-11
Substrate A	101159GR	QNSA220620	2024-JUN-20
Substrate B+	101120GR	QNSB+220620	2024-JUN-20
Wash Buffer	101158GR	QNSW220615	2024-JUN-15

Spot #	Analyte	Units	Cal Conc.	ULOQ	LLOQ	LLD
1	IFN γ	pg/mL	7130	N/A	N/A	5.40
2	IL-1 β	pg/mL	12950	N/A	N/A	12.41
3	IL-6	pg/mL	4800	N/A	N/A	2.90
4	TNF α	pg/mL	2980	N/A	N/A	3.40
5	REF	N/A	N/A	N/A	N/A	N/A

Quality Release:

for Research Use Only

101896GR 2.0