


Q-Plex™ ARRAY

CERTIFICATE OF ANALYSIS

Product Name: Q-Plex™ Human
 Custom (3-Plex)
 Order No. 20977
 REF Part No. 107749GR
 LOT Lot No. HCUK220531
 Exp Date: 2024-JAN-07
 Software HCUMA220531-UW
 Product Code:



Reconstitute Calibrator with 650 µL of Complete Sample Diluent

Component	Part No.	Lot No.	Exp. Date
Plate	101852GR	HCUP220906	2024-SEP-06
Sample Diluent	101104HU	HQ6S220107	2024-JAN-07
Calibrator	101855GR	HCUM220531	2024-SEP-14
Competitor	101890GR	HCUC191111	2024-SEP-14
Detection	101856GR	HCUD220909	2024-SEP-09
Streptavidin HRP	101173GR	QNSN220715	2024-JAN-15
Substrate A	101159GR	QNSA220713	2024-JUL-13
Substrate B+	101120GR	QNSB+220713	2024-JUL-13
Wash Buffer	101158GR	QNSW220803	2024-AUG-03

Spot #	Analyte	Units	Cal	ULOO	LLOO	LLD
1	CRP	mg/mL	27.18	27.18	0.11	0.055
2	IL-6	pg/mL	346.15	346.15	0.47	0.24
3	TNFa	pg/mL	719.87	719.87	0.59	0.30
4	REF	N/A	N/A	N/A	N/A	N/A

Quality Release:

[Signature]
 for Research Use Only

101896GR 2.0