


# Q-Plex™ ARRAY

## CERTIFICATE OF ANALYSIS

Product Name: Q-Plex™ Mouse  
 Custom (3-Plex)  
 Order No. 21669D  
 Part No. 107749GR  
 Lot No. MCUK191126 18395 001  
 Exp Date: 2023-NOV-20  
 Software MA1M191126-CC3  
 Product Code: MA1M191126-CC3



**Reconstitute Calibrator with 550 µL of Sample Diluent**

| Component        | Part No. | Lot No.     | Exp. Date   |
|------------------|----------|-------------|-------------|
| Plate            | 101852GR | MCUP220609  | 2024-JUN-09 |
| Sample Diluent   | 110957MS | MCSS220607  | 2024-JUN-07 |
| Calibrator       | 102655MS | MA1M191126  | 2024-SEP-12 |
| Competitor       | N/A      | N/A         | N/A         |
| Detection        | 101856GR | MCUD220610  | 2024-JUN-10 |
| Streptavidin HRP | 101173GR | QNSN220520  | 2023-NOV-20 |
| Substrate A      | 101159GR | QNSA211217  | 2023-DEC-17 |
| Substrate B+     | 101120GR | QNSB+211217 | 2023-DEC-17 |
| Wash Buffer      | 101158GR | QNSW220422  | 2024-APR-22 |

| Spot # | Analyte | Units | Cal   | ULOO | LLOO | LLD   |
|--------|---------|-------|-------|------|------|-------|
| 1      | IL-1b   | pg/mL | 12950 | N/A  | N/A  | 12.41 |
| 2      | IL-6    | pg/mL | 4800  | N/A  | N/A  | 2.90  |
| 3      | TNFa    | pg/mL | 2980  | N/A  | N/A  | 3.40  |
| 4      | REF     | N/A   | N/A   | N/A  | N/A  | N/A   |

Quality Release:

*[Handwritten Signature]*

101896GR 2.0

**for Research Use Only**