


## CERTIFICATE OF ANALYSIS

Product Name: Q-Plex™ Mouse  
 Custom (1-Plex)  
 Order No. 21669D  
**REF** Part No. 107749GR  
**LOT** Lot No. MCUK191126 18395 002  
 Exp Date: 2024-JAN-15  
**Software Product Code: MA1M191126-MCP1**



### Reconstitute Calibrator with 550 µL of Sample Diluent

Component	<b>REF</b> Part No.	<b>LOT</b> Lot No.	Exp. Date
Plate	101852GR	MCUP220808	2024-AUG-08
Sample Diluent	110957MS	MCSS220602	2024-MAY-01
Calibrator	102655MS	MA1M191126	2024-SEP-12
Competitor	N/A	N/A	N/A
Detection	101856GR	MCUD220810	2024-AUG-10
Streptavidin HRP	101173GR	QNSN220715	2024-JAN-15
Substrate A	101159GR	QNSA220713	2024-JUL-13
Substrate B+	101120GR	QNSB+220713	2024-JUL-13
Wash Buffer	101158GR	QNSW220803	2024-AUG-03

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	MCP-1	pg/mL	3150	N/A	N/A	3.40
2	REF	N/A	N/A	N/A	N/A	N/A

Quality Release:

101896GR 2.0

for Research Use Only