


# Q-Plex™ ARRAY

## CERTIFICATE OF ANALYSIS

Product Name: Q-Plex™ Human  
MMP (6-Plex)  
Order No. 21792  
Part No. 340949HU  
Lot No. MMPK220518 18381 001  
Exp Date: 2023-JUL-16  
Software  
Product Code: **MMPM220518**



**Reconstitute Calibrator with 300 µL of Sample Diluent**

Component	Part No.	Lot No.	Exp. Date
Plate	340952HU	MMPP220630	2024-JUN-30
Sample Diluent	340957HU	MMPS220124	2024-JAN-24
Calibrator	340955HU	MMPM220518	2024-SEP-09
Assay Diluent	340983HU	MMPV220124	2024-JAN-24
Detection	340956HU	MMPD210716	2023-JUL-16
Streptavidin HRP	101173GR	QNSN220715	2024-JAN-15
Substrate A	101159GR	QNSA220713	2024-JUL-13
Substrate B+	101120GR	QNSB+220713	2024-JUL-13
Wash Buffer	101158GR	QNSW220803	2024-AUG-03

Spot #	Analyte	Units	Cal	ULOQ	LLOQ	LLD
1	MMP 1	pg/mL	12563	12563	51.70	25.85
2	MMP 2	pg/mL	85197	85197	350.60	175.30
3	MMP 3	pg/mL	4579	4579	18.84	9.42
4	MMP 7	pg/mL	1745	1745	7.18	3.59
5	MMP 9	pg/mL	44914	44914	184.83	92.42
6	MMP 13	pg/mL	32943	32943	135.57	67.78
7	REF	N/A	N/A	N/A	N/A	N/A

Quality Release:

  
*for Research Use Only*

101896GR 2.0