


CERTIFICATE OF ANALYSIS

Product Name: Q-Plex™ Porcine Cytokine - High Sensitivity (4-Plex)
 Order No. 22163
 Part No. 119149PC
 Lot No. PCSK170306 18418 001
 Exp Date: 2023-JUN-15
 Software Product Code: **PCSM170306CY**



Reconstitute Calibrator with 500 µL of Sample Diluent

Component	REF	Part No.	LOT	Lot No.	Exp. Date
Plate		119152PC		PCSP210726	2023-JUL-26
Sample Diluent		119157PC		PCSS220321	2024-MAR-21
Calibrator		119155PC		PCSM170306	2024-SEP-16
Competitor		N/A		N/A	N/A
Detection		119156PC		PCSD210615	2023-JUN-15
Streptavidin HRP		101173GR		QNSN220715	2024-JAN-15
Substrate A		101159GR		QNSA220713	2024-JUL-13
Substrate B+		101120GR		QNSB+220713	2024-JUL-13
Wash Buffer		101158GR		QNSW220803	2024-AUG-03

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	IL-1b	pg/mL	1600	N/A	N/A	2.10
2	IL-6	pg/mL	2000	N/A	N/A	2.50
3	IL-8	pg/mL	2300	N/A	N/A	2.90
4	TNFa	pg/mL	1350	N/A	N/A	1.90
5	REF	N/A	N/A	N/A	N/A	N/A

Quality Release:

Handwritten signature: Nancy Fisher

101896GR 2.0

for Research Use Only