


CERTIFICATE OF ANALYSIS

 Q-Plex™ Human Malaria (5-Plex)
 Software Product Code: HMRM191125-ML2

☐ **REF** Part Number: 565949HU
☐ **LOT** Lot Number: HMRK191125 18498 001
☒ Expiration Date: 2023-MAY-15

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	565952HU	HMRP220801	N/A	2024-FEB-01
Sample Diluent	N/A	N/A	N/A	N/A
Calibrator Diluent	565928HU	HMRU211115	6 mL	2023-MAY-15
Assay Diluent	565983HU	HMRV211116	N/A	2023-MAY-16
Competitor	565990HU	HMRC191014	6 mL	2024-APR-24
Calibrator	565955HU	HMRM191125	350 µL	2024-APR-24
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	565956HU	HMRD211119	N/A	2023-MAY-19
Streptavidin HRP	101173GR	QNSN220815	N/A	2024-FEB-15
Substrate A+	105959GR	MOSA+220810	N/A	2024-AUG-10
Substrate B	105920GR	MOSB220810	N/A	2024-AUG-10
Wash Buffer	101158GR	QNSW220909	N/A	2024-SEP-09

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	HRP2	pg/mL	1890	2800	0.68	0.53
2	LDH-Pan	pg/mL	72889	67000	16.36	7.28
3	LDH-Pv	pg/mL	18063	19200	4.69	4.23
4	LDH-Pf	pg/mL	21418	20800	5.08	2.33
5	*CRP	ng/mL	45286	38000	9.28	7.65
6	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS
 *Competitive Assay


QC Release:



Date:

24OCT22

QA Release:



Date:

24OCT22

for Research Use Only