

CERTIFICATE OF ANALYSIS

 Q-Plex™ Human Pneumococcal IgG (9-Plex)
 Software Product Code: HPNM007SPWHO9

LOT
REF

Part Number: 481649HU


Lot Number: HPNK220926 18438 001



Expiration Date: 2023-AUG-16

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	481652HU	HPNP220816	N/A	2023-AUG-16
Sample Diluent	481157HU	HPNS220926	N/A	2023-SEP-26
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	N/A	N/A	N/A	N/A
Calibrator	481155HU	HPNM220926	500 µL	2023-OCT-06
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	481156HU	HPND220823	N/A	2023-AUG-23
Streptavidin HRP	N/A	N/A	N/A	N/A
Substrate A	101159GR	QNSA220825	N/A	2024-AUG-25
Substrate B+	101120GR	QNSB+220825	N/A	2024-AUG-25
Wash Buffer	101158GR	QNSW220909	N/A	2024-SEP-09

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	2	µg/mL	2.46	0.434	0.00078	0.00039
2	10A	µg/mL	1.30	0.0567	0.000416	0.000208
3	11A	µg/mL	0.51	0.091	0.000156	0.000078
4	15B	µg/mL	1.69	0.294	0.000533	0.000267
5	17F	µg/mL	0.85	0.140	0.00108	0.00054
6	19A	µg/mL	1.39	0.245	0.00182	0.00091
7	20	µg/mL	1.05	0.182	0.000338	0.000169
8	22F	µg/mL	0.95	0.168	0.001209	0.000605
9	33F	µg/mL	1.07	0.0469	0.000338	0.000169
10	CWPS	N/A	N/A	N/A	N/A	N/A
11	DK anti-Hu	N/A	N/A	N/A	N/A	N/A


Q-Plex™ ARRAYS

QC Release:

Date: 06 OCT 22

QA Release:

Date: 06 OCT 22

for Research Use Only