

101896GR - 4.2



Quansys Biosciences
365 N 600 W Logan, UT 84321
(435) 752-0531


CERTIFICATE OF ANALYSIS

Q-Plex™ Human Custom 8-Plex
Software Product Code: HA1M200518-NJH

REF Part Number: 107749GR
LOT Lot Number: HCUK200518 18646 001
Expiration Date: 2024-JAN-07

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	101852GR	HCUP221110	N/A	2024-NOV-10
Sample Diluent	101104HU	HQ6S220107	N/A	2024-JAN-07
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	N/A	N/A	N/A	N/A
Calibrator	102655HU	HA1M200518	650 µL	2024-NOV-11
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	110956HU 01	HCSDQ6220407	N/A	2024-APR-07
Streptavidin	101173GR	QNSN221024	N/A	2024-APR-24
Substrate A	101159GR	QNSA221020	N/A	2024-OCT-20
Substrate B+	101120GR	QNSB+221020	N/A	2024-OCT-20
Wash Buffer	101158GR	QNSW220909	N/A	2024-SEP-09

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	IFN γ	pg/mL	8000	8000	10.97	5.49
2	IL-2	pg/mL	4600	4600	6.31	3.16
3	IL-4	pg/mL	1600	1600	2.19	1.10
4	IL-5	pg/mL	3100	3100	4.25	2.13
5	IL-6	pg/mL	3100	3100	4.25	2.13
6	IL-10	pg/mL	6800	6800	9.33	4.66
7	IL-12p70	pg/mL	2800	2800	3.84	1.92
8	TNF α	pg/mL	5400	5400	7.41	3.70
9	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: [Signature] Date: 15 NOV 22
 QA Release: [Signature] Date: 15 NOV 22

for Research Use Only