

QUANSYS
 BIOSCIENCES
CERTIFICATE OF ANALYSIS


Quansys Biosciences
 365 N 600 W Logan, UT 84321
 (435) 752-0531

Q-Plex™ NHP Custom 13-Plex
 Software Product Code: NHPM220606-MD13

REF Part Number: 107749GR
LOT Lot Number: NHPCUK220606 18592 001
 Expiration Date: 2023-OCT-18

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	101852GR	NHPCUP221018	N/A	2023-OCT-18
Sample Diluent	101104HU	HQ6S220909	N/A	2024-SEP-09
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	N/A	N/A	N/A	N/A
Calibrator	109255HU	NHPM220606	500 µL	2023-NOV-01
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	101856GR	NHPCUD221021	6 mL	2023-OCT-21
Streptavidin HRP	101173GR	QNSN220715	N/A	2024-JAN-15
Substrate A	101159GR	QNSA220825	N/A	2024-AUG-25
Substrate B+	101120GR	QNSB+220825	N/A	2024-AUG-25
Wash Buffer	101158GR	QNSW220909	N/A	2024-SEP-09

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	GMCSF	pg/mL	5970	5970	8.19	4.09
2	GROa	pg/mL	1760	1760	2.41	1.21
3	IFN γ	pg/mL	820	820	1.12	0.56
4	IL-1b	pg/mL	3660	3660	5.02	2.51
5	IL-6	pg/mL	910	910	1.25	0.62
6	IL-8	pg/mL	1920	1920	2.63	1.32
7	IL-10	pg/mL	2460	2460	3.37	1.69
8	IL-12p70	pg/mL	2630	2630	3.61	1.80
9	IP-10	pg/mL	2910	970	3.99	2.00
10	MCP-1	pg/mL	1880	1880	7.74	3.87
11	MIG	pg/mL	12370	12370	16.97	8.48
12	MIP-1a	pg/mL	1800	600	2.47	1.23
13	TNFA	pg/mL	190	190	0.26	0.13
14	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: Sophia Adams

Date: 01 NOV 22

QA Release: R. Wa

Date: 01 NOV 22

for Research Use Only