


CERTIFICATE OF ANALYSIS

 Q-Plex™ Human Angiogenesis (9-Plex)
 Software Product Code: HA1M200518-AGQ


REF Part Number: 150233HU
LOT Lot Number: HAGK200518 18711 001
EXP Expiration Date: 2023-OCT-13

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	150252HU	HAGP220815	N/A	2024-AUG-15
Sample Diluent	101104HU	HQ6S220107	N/A	2024-JAN-07
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	N/A	N/A	N/A	N/A
Calibrator	102655HU	HA1M200518	650 µL	2024-DEC-01
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	150256HU_01	HAGDQ6220302	N/A	2024-MAR-02
Streptavidin HRP	101173GR	QNSN221024	N/A	2024-APR-24
Substrate A	101159GR	QNSA221020	N/A	2024-OCT-20
Substrate B+	101120GR	QNSB+221020	N/A	2024-OCT-20
Wash Buffer	101158GR	QNSW211013	N/A	2023-OCT-13

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	ANG-2	pg/mL	8600	8600	11.80	5.90
2	FGF	pg/mL	5200	5200	7.13	3.57
3	HGF	pg/mL	13700	13700	18.79	9.40
4	IL-8	pg/mL	1800	1800	2.47	1.23
5	PDGF	pg/mL	4300	4300	5.90	2.95
6	TIMP-1	pg/mL	14200	14200	19.48	9.74
7	TIMP-2	pg/mL	18200	18200	24.97	12.48
8	TNFα	pg/mL	5400	5400	7.41	3.70
9	VEGF	pg/mL	2000	2000	2.74	1.37
10	REF	N/A	N/A	N/A	N/A	N/A


Q-Plex™ ARRAYS


QC Release:



Date:

02 DEC 22

QA Release:



Date:

02 DEC 22

for Research Use Only