


CERTIFICATE OF ANALYSIS

 Q-Plex™ Human Custom 1-Plex
 Software Product Code: HA1M200518-PDGF

LOT
REF Part Number: 107749GR
 Lot Number: HCUK200518 18759 001
 Expiration Date: 2024-MAR-02

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	101852GR	HCUP221130	N/A	2024-NOV-30
Sample Diluent	101104HU	HQ6S220909	N/A	2024-SEP-09
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	N/A	N/A	N/A	N/A
Calibrator	102655HU	HA1M200518	650 µL	2024-DEC-13
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	150256HU 01	HAGDQ6220302	N/A	2024-MAR-02
Streptavidin HRP	101173GR	QNSN221024	N/A	2024-APR-24
Substrate A	101159GR	QNSA221020	N/A	2024-OCT-20
Substrate B+	101120GR	QNSB+221020	N/A	2024-OCT-20
Wash Buffer	101158GR	QNSW221123	N/A	2024-NOV-23

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	PDGF-BB	pg/mL	4300	4300	5.90	2.95
2	PDGF-BB	pg/mL	4300	4300	5.90	2.95
3	REF	N/A	N/A	N/A	N/A	N/A


Q-Plex™ ARRAYS

QC Release:

Sophia Adams

Date:

13 DEC 22

QA Release:

[Signature]

Date:

13 DEC 22

for Research Use Only