

101896GR - 5.0



Quansys Biosciences
365 N 600 W Logan, UT 84321
(435) 752-0531


CERTIFICATE OF ANALYSIS

Q-Plex™ Human Angiogenesis (9-Plex)
Software Product Code: HA1M200518-AG

REF Part Number: 150233HU
LOT Lot Number: HAGK200518 19003 001
Expiration Date: 2024-JUN-28

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	150252HU	HAGP220815	N/A	2024-AUG-15
Sample Diluent	101104HU	HQ6S220909	N/A	2024-SEP-09
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	N/A	N/A	N/A	N/A
Calibrator	102655HU	HA1M200518	650 µL	2025-FEB-14
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	150256HU_01	HAGDQ6221230	N/A	2024-DEC-30
Streptavidin HRP	101173GR	QNSN221228	N/A	2024-JUN-28
Substrate A	101159GR	QNSA221221	N/A	2024-DEC-21
Substrate B+	101120GR	QNSB+221221	N/A	2024-DEC-21
Wash Buffer	101158GR	QNSW221228	N/A	2024-DEC-28

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	ANG-2	pg/mL	8600	8600	11.80	5.90
2	FGF	pg/mL	5200	1733	7.13	3.57
3	HGF	pg/mL	13700	13700	18.79	9.40
4	IL-8	pg/mL	1800	1800	2.47	1.23
5	PDGF	pg/mL	4300	4300	5.90	2.95
6	TIMP-1	pg/mL	14200	14200	58.44	29.22
7	TIMP-2	pg/mL	18200	18200	24.97	12.48
8	TNFα	pg/mL	5400	5400	7.41	3.70
9	VEGF	pg/mL	2000	2000	2.74	1.37
10	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: *Erin Adams* Date: 14 FEB 23
 QA Release: *J. Wa* Date: 14 FEB 23
 for Research Use Only