



NHP Multiplex ELISA Kits

The four Q-Plex™ NHP kits encompass over 35 cytokines, chemokines, and cancer biomarkers. Each Q-Plex NHP kit has a customizable list of biomarkers that can be used for many types of research projects. Choose the biomarkers most suitable to your research and get custom plates with pre-printed capture material made just for you.

The Q-Plex NHP ELISA kits provides researchers with a fast and easy-to-use means of generating biomarker profiles for cytokines, chemokines, and cancer biomarkers for up to 80 samples on one plate!

Catalog #	Product	Assays				
359233MB	Q-Plex™ NHP Anti-Inflammatory Panel 1 (7-Plex)	IL-1Rα IL-2Rα	IL-4 IL-6R	IL-10 IL-13	TNFRII	
189233MB	Q-Plex™ NHP Cancer Panel 1 (6-Plex)	Ang-2 CRP	PDGF TIMP1	TIMP2 VEGF		
120933MB	Q-Plex™ NHP Chemokine Screen (9-Plex)	CXCL-1 IL-8	IP-10 MCP-1	MIF MIG	MIP-1α MIP-1β	RANTES
110933MB	Q-Plex™ NHP Cytokine Screen (16-Plex)	GCSF GMCSF IFNγ IL-1β	IL-2 IL-4 IL-5 IL-6	IL-8 IL-10 IL-12p70 IL-13	IL-15 IL-17a IL-18 TNFα	

Request an Analyte

If you don't see the analyte you need, request it! At Quansys, we want to make sure you have what you need when it comes to our Q-Plex kits. Requesting an analyte can be done over the phone, via email, or through an online submission form. Be sure to provide any information about the analytes you want along with any other details about your project or samples tested with it.

(435) 752-0531 | info@quansybio.com | <https://www.quansysbio.com/request-an-analyte/>

Visit <https://www.quansysbio.com/> to learn more about our available assays, options for building your own assay, custom development, and more.



Quansys Biosciences

365 North 600 West
Logan, UT 84321
www.quansysbio.com

Fax: 435-750-6869
Phone: 435-752-0531
Toll Free: 888-QUANSYS (782-6797)



Quansys is an ISO 9001:2015 and ISO 13485:2016 registered company and complies with GMP. Products are designed, developed, and manufactured according to the procedures outlined in our Quality Management System.

©2022 Quansys Biosciences Inc. All rights reserved. Q-View and Q-Plex are Trademarks of Quansys Biosciences Inc. ISO is a trademark of the International Standards Organization. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.