

101896GR - 6.0



Quansys Biosciences  
365 N 600 W Logan, UT 84321  
(435) 752-0531

### CERTIFICATE OF ANALYSIS


Q-Plex™ Human Custom 3-Plex  
Software Product Code: HA2M200303-IFN

**REF** Part Number: 107749GR  
**LOT** Lot Number: HCUK220712 19465 001  
Expiration Date: 2024-AUG-16

**Reconstitute Calibrator 1 using 600 µL of Sample Diluent**  
**Reconstitute Calibrator 2 using 400 µL of Calibrator 1**

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	101852GR	HCUP230608	N/A	2025-JUN-08
Sample	101104HU	HQ6S230525	N/A	2025-MAY-25
Calibrator	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Calibrator 1	102655HU	HA1M220712	600	2025-JUN-19
Calibrator 2	102855HU	HA2M200303	400	2025-JUN-19
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	101856GR	HCUD230609	6 mL	2025-JUN-09
Streptavidin	101173GR	QNSN230216	N/A	2024-AUG-16
Substrate A	101159GR	QNSA230302	N/A	2025-MAR-02
Substrate B+	101120GR	QNSB+230302	N/A	2025-MAR-02
Wash Buffer	101158GR	QNSW230314	N/A	2025-MAR-14

Spot #	Analyte	Units	Cal Conc	ULOQ	LLOQ	LLD
1	IFN $\alpha$	pg/mL	385	385	0.53	0.26
2	IFN $\beta$	pg/mL	21700	21700	29.77	14.88
3	IFN $\gamma$	pg/mL	3990	3990	5.47	2.74
4	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: [Signature]  
QA Release: [Signature]  
for Research Use Only

Date: 19 JUN 23  
Date: 19 JUN 23