

101896GR - 7.0

QUANSYS
BIOSCIENCES
CERTIFICATE OF ANALYSIS

Quansys Biosciences
365 N 600 W Logan, UT 84321
(435) 752-0531


Q-Plex™ Human Pneumococcal IgG (9-Plex)
Software Product Code: HPNP007SPWHO9



REF Part Number: 481649HU
LOT Lot Number: HPNK220926 19560 001
Expiration Date: 2024-JUN-12

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	481652HU	HPNP230612	N/A	2024-JUN-12
Sample Diluent	481157HU	HPNS230621	N/A	2024-JUN-21
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	N/A	N/A	N/A	N/A
Calibrator	481155HU	HPNM220926	500 µL	2024-JUL-06
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	481156HU	HPND230622	N/A	2024-JUN-22
Streptavidin HRP	N/A	N/A	N/A	N/A
Substrate A	101159GR	QNSA230302	N/A	2025-MAR-02
Substrate B+	101120GR	QNSB+230302	N/A	2025-MAR-02
Wash Buffer	101158GR	QNSW230314	N/A	2025-MAR-14

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	2	2.46	µg/mL	0.43	0.00078	0.00039
2	10A	1.30	µg/mL	0.057	0.00042	0.00021
3	11A	0.51	µg/mL	0.091	0.00016	0.000078
4	15B	1.69	µg/mL	0.29	0.00053	0.00027
5	17F	0.85	µg/mL	0.14	0.0011	0.00054
6	19A	1.39	µg/mL	0.25	0.0018	0.00091
7	20	1.05	µg/mL	0.18	0.00034	0.00017
8	22F	0.95	µg/mL	0.17	0.0012	0.00061
9	33F	1.07	µg/mL	0.047	0.00034	0.00017
10	CWPS	N/A	N/A	N/A	N/A	N/A
11	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: [Signature]
QA Release: [Signature]
for Research Use Only

Date: 11 JUL 23
Date: 11 JUL 23