

101896GR - 7.0



Quansys Biosciences
365 N 600 W Logan, UT 84321
(435) 752-0531

CERTIFICATE OF ANALYSIS

Q-Plex™ Human Custom Stripwell 14-Plex
Software Product Code: HA1M220712-GS14




REF Part Number: 107949GR
LOT Lot Number: HCUK220712 19458 001
Expiration Date: 2024-DEC-28

Reconstitute Calibrator 1 with 600 µL of Sample Diluent
Reconstitute Calibrator 2 with 400 µL of Calibrator 1
Reconstitute Calibrator 3 with 300 µL of Calibrator 2

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	101852GR	HCUP230725	N/A	2025-JUL-25
Sample Diluent	101104HU	HQ6S230525	N/A	2025-MAY-25
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Calibrator 1	102655HU	HA1M220712	600 µL	2025-AUG-28
Calibrator 2	102855HU	HA2M200303	400 µL	2025-AUG-28
Calibrator 3	101855GR	HCUM230705	300 µL	2025-AUG-28
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	101856GR	HCUD230731	N/A	2025-JUL-31
Streptavidin HRP	101173GR	QNSN230628	N/A	2024-DEC-28
Substrate A	101159GR	QNSA230602	N/A	2025-JUN-02
Substrate B+	101120GR	QNSB+230602	N/A	2025-JUN-02
Wash Buffer	101158GR	QNSW230718	N/A	2025-JUL-18

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	EGF	319	pg/mL	319	0.44	0.22
2	FGF	4950	pg/mL	4950	6.79	3.40
3	GCSF	4420	pg/mL	4420	6.06	3.03
4	GMCSF	6010	pg/mL	6010	8.24	4.12
5	HGF	15160	pg/mL	15160	20.80	10.40
6	Leptin	57	ng/mL	19	0.23	0.12
7	MCP-1	1720	pg/mL	573	7.08	3.54
8	MCP-2	585	pg/mL	585	7.22	3.61
9	MCP-3	114	pg/mL	114	0.16	0.078
10	PAI-1	2675	pg/mL	2675	3.67	1.83
11	Resistin	11100	pg/mL	11100	15.23	7.61
12	TIMP-1	12840	pg/mL	12840	52.84	26.42
13	TIMP-2	19135	pg/mL	19135	26.25	13.12
14	VEGF	1890	pg/mL	1890	2.59	1.30
15	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: [Signature] Date: 28 AUG 23

QA Release: [Signature] Date: 28 Aug 23

for Research Use Only