

101896GR - 7.0

QUANSYS
BIOSCIENCES
CERTIFICATE OF ANALYSIS

Quansys Biosciences
365 N 600 W Logan, UT 84321
(435) 752-0531


Q-Plex™ Human Cytokine Inflammation Panel 1 (9-Plex)
Software Product Code: HA1M220712-IF



REF Part Number: 110433HU
LOT Lot Number: HClK220712 20158 001
Expiration Date: 2025-MAR-02

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	110452HU	HCIP230606	N/A	2025-JUN-06
Sample Diluent	101104HU	HQ6S230525	N/A	2025-MAY-25
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	N/A	N/A	N/A	N/A
Calibrator	102655HU	HA1M220712	600 µL	2025-NOV-30
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	110956HU_01	HCSDQ6230612	N/A	2025-JUN-12
Streptavidin HRP	101173GR	QNSN230906	N/A	2025-MAR-06
Substrate A	101159GR	QNSA230302	N/A	2025-MAR-02
Substrate B+	101120GR	QNSB+230302	N/A	2025-MAR-02
Wash Buffer	101158GR	QNSW230919	N/A	2025-SEP-19

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	IL-1 α	3940	pg/mL	3940	5.40	2.70
2	IL-1 β	11800	pg/mL	11800	16.19	8.09
3	IL-2	5150	pg/mL	5150	7.06	3.53
4	IL-4	1670	pg/mL	1670	2.29	1.15
5	IL-6	2850	pg/mL	2850	3.91	1.95
6	IL-8	2080	pg/mL	2080	2.85	1.43
7	IL-10	8400	pg/mL	8400	11.52	5.76
8	IFN γ	9775	pg/mL	9775	13.41	6.70
9	TNF α	6020	pg/mL	6020	8.26	4.13
10	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release:

Date:

30 NOV 23

QA Release:

Date:

30 NOV 23

for Research Use Only