

101896GR - 7.0

**QUANSYS**  
BIOSCIENCES  
**CERTIFICATE OF ANALYSIS**

Quansys Biosciences  
365 N 600 W Logan, UT 84321  
(435) 752-0531


Q-Plex™ Human MMP (6-Plex)  
Software Product Code: MMPM220518



**REF** Part Number: 340949HU  
**LOT** Lot Number: MMPK220518 20218 001  
Expiration Date: 2024-OCT-13

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	340952HU	MMPP230828	N/A	2025-AUG-28
Sample Diluent	340957HU	MMPS230824	N/A	2025-AUG-24
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	340983HU	MMPV230824	N/A	2025-AUG-24
Competitor	N/A	N/A	N/A	N/A
Calibrator	340955HU	MMPM220518	300 µL	2025-DEC-14
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	340956HU	MMPD230926	N/A	2025-SEP-26
Streptavidin HRP	101173GR	QNSN231025	N/A	2025-APR-25
Substrate A	101159GR	QNSA231205	N/A	2024-OCT-13
Substrate B+	101120GR	QNSB+231205	N/A	2024-OCT-13
Wash Buffer	101158GR	QNSW230718	N/A	2025-JUL-18

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	MMP 1	12563	pg/mL	12563.00	51.70	25.85
2	MMP 2	85197	pg/mL	85197.00	350.60	175.30
3	MMP 3	4579	pg/mL	4579.00	18.84	9.42
4	MMP 7	1745	pg/mL	1745.00	7.18	3.59
5	MMP 9	44914	pg/mL	44914.00	184.83	92.42
6	MMP 13	32943	pg/mL	32943.00	135.57	67.78
7	REF	N/A	N/A	N/A	N/A	N/A



**Q-Plex™** ARRAYS

QC Release: [Signature]  
QA Release: [Signature]  
for Research Use Only

Date: 20 DEC 23  
Date: 20 DEC 23