

101896GR - 7.0

# QUANSYS BIOSCIENCES CERTIFICATE OF ANALYSIS

Quansys Biosciences  
365 N 600 W Logan, UT 84321  
(435) 752-0531


Q-Plex™ Human Female Hormone (6-Plex)  
Software Product Code: FHCM190715-6P



**REF** Part Number: 332149HU  
**LOT** Lot Number: FHCK190715 20154 001  
Expiration Date: 2024-OCT-13

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	332152HU	FHCP231201	N/A	2024-DEC-01
Sample Diluent	101122HU	HQ16S231221	N/A	2024-DEC-21
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	332190HU	FHCC231208	N/A	2024-DEC-08
Calibrator	332155HU	FHCM190715	550 µL	2025-JUN-21
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	332156HU	FHCD231207	N/A	2024-DEC-07
Streptavidin HRP	101173GR	QNSN231025	N/A	2025-APR-25
Substrate A	101159GR	QNSA231205	N/A	2024-OCT-13
Substrate B+	101120GR	QNSB+231205	N/A	2024-OCT-13
Wash Buffer	101158GR	QNSW231024	N/A	2025-OCT-24

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	Adiponectin	25	ng/mL	25	0.10	0.051
2	*Cortisol	714	ng/mL	714	0.98	0.49
3	C-Peptide	36	ng/mL	36	0.049	0.025
4	*E1G	522	ng/mL	522	0.72	0.36
5	FSH	193	mIU/mL	193	0.79	0.40
6	hCGβ	40	ng/mL	40	0.055	0.027
7	REF	N/A	N/A	N/A	N/A	N/A



**Q-Plex™ ARRAYS**  
\*Competitive Assay

QC Release:

*[Signature]*

Date:

21 DEC 23

QA Release:

*[Signature]*

Date:

21 DEC 23

for Research Use Only