

101896GR - 7.0

QUANSYS BIOSCIENCES CERTIFICATE OF ANALYSIS

Quansys Biosciences
365 N 600 W Logan, UT 84321
(435) 752-0531

Q-Plex™ Human Type 2 Immunity (5-Plex)
Software Product Code: HA1M220712-T2




REF Part Number: 114233HU
LOT Lot Number: H2IK220712 20253 001
Expiration Date: 2024-OCT-13

Reconstitute Calibrator 1 using 600 µL of Sample Diluent
Reconstitute Calibrator 2 using 500 µL of Calibrator 1

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	114252HU	H2IP231214	N/A	2025-DEC-14
Sample Diluent	101104HU	HQ6S230525	N/A	2025-MAY-25
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Calibrator 1	102655HU	HA1M220712	600 µL	2025-DEC-20
Calibrator 2	118655HU	HCFM171206	500 µL	2025-DEC-20
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	101856GR	HCUD231219	N/A	2025-DEC-19
Streptavidin HRP	101173GR	QNSN231025	N/A	2025-APR-25
Substrate A	101159GR	QNSA231205	N/A	2024-OCT-13
Substrate B+	101120GR	QNSB+231205	N/A	2024-OCT-13
Wash Buffer	101158GR	QNSW231122	N/A	2025-NOV-22

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	IL-4	1670	pg/mL	1670	2.29	1.15
2	IL-5	2950	pg/mL	2950	4.05	2.02
3	IL-10	8400	pg/mL	8400	11.52	5.76
4	IL-13	2865	pg/mL	2865	3.93	1.97
5	IL-33	315	pg/mL	315	0.43	0.22
6	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: [Signature] Date: 03 JAN 24
 QA Release: [Signature] Date: 03 JAN 24
 for Research Use Only