

101896GR - 7.0

QUANSYS
BIOSCIENCES
CERTIFICATE OF ANALYSIS

Quansys Biosciences
365 N 600 W Logan, UT 84321
(435) 752-0531


Q-Plex™ Human Innate Immunity (7-Plex)
Software Product Code: HA1M200518-IN

REF Part Number: 111233HU
LOT Lot Number: HINK200518 20410 001
Expiration Date: 2024-DEC-01

Reconstitute Calibrator 1 using 650 µL of Sample Diluent
Reconstitute Calibrator 2 using 400 µL of Calibrator 1

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	111252HU	HINP221201	N/A	2024-DEC-01
Sample Diluent	101104HU	HQ6S230525	N/A	2025-MAY-25
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Calibrator 1	102655HU	HA1M200518	650 µL	2025-MAY-31
Calibrator 2	102855HU	HA2M200303	400 µL	2026-JAN-31
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	101856GR	HCUD221209	6 mL	2024-DEC-09
Streptavidin HRP	101173GR	QNSN230906	N/A	2025-MAR-06
Substrate A	101159GR	QNSA231006	N/A	2025-MAR-13
Substrate B+	101120GR	QNSB+231006	N/A	2025-MAR-13
Wash Buffer	101158GR	QNSW231218	N/A	2025-DEC-18

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	IL-1α	3900	pg/mL	3900	5.35	2.67
2	IL-1β	11800	pg/mL	11800	16.19	8.09
3	IFNα	380	pg/mL	380	0.52	0.26
4	IFNβ	21700	pg/mL	21700	29.77	14.88
5	IFNγ	8000	pg/mL	8000	10.97	5.49
6	GMCSF	6000	pg/mL	6000	8.23	4.12
7	IL-18	235	pg/mL	235	0.32	0.16
8	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: [Signature]

Date: 31 JAN 24

QA Release: [Signature]

Date: 31 JAN 24

for Research Use Only