

101896GR - 7.0

QUANSYS BIOSCIENCES

CERTIFICATE OF ANALYSIS

Quansys Biosciences
365 N 600 W Logan, UT 84321
(435) 752-0531


Q-Plex™ Mouse Cytokine Panel 1 HS (6-Plex)
Software Product Code: MHSM230927-P16P



REF Part Number: 111349MS
LOT Lot Number: MCP1K230927 20434 001
Expiration Date: 2024-JUN-30

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	111352MS	MCP1P240212	N/A	2025-FEB-12
Sample Diluent	112457MS	MHSS230512	N/A	2024-JUN-30
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	N/A	N/A	N/A	N/A
Calibrator	112455MS	MHSM230927	550 µL	2025-FEB-14
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	112456MS	MHSD231013	6 mL	2024-OCT-13
Streptavidin HRP	101173GR	QNSN230216	N/A	2024-AUG-16
Substrate A	101159GR	QNSA231006	N/A	2025-MAR-13
Substrate B+	101120GR	QNSB+231006	N/A	2025-MAR-13
Wash Buffer	101158GR	QNSW240103	N/A	2026-JAN-03

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	IFN γ	800	pg/mL	800	1.10	0.55
2	IL-1 α	620	pg/mL	620	0.85	0.43
3	IL-1 β	2390	pg/mL	2390	3.28	1.64
4	IL-6	1110	pg/mL	1110	1.52	0.76
5	IL-10	520	pg/mL	520	0.71	0.36
6	TNF α	390	pg/mL	390	0.53	0.27
7	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release:

[Signature]

Date:

14 FEB 24

QA Release:

[Signature]

Date:

14 FEB 24

for Research Use Only