

101896GR - 7.0



Quansys Biosciences
365 N 600 W Logan, UT 84321
(435) 752-0531

CERTIFICATE OF ANALYSIS


Q-Plex™ NHP Custom Stripwell 18-Plex
Software Product Code: NHPM220606-TB18

REF Part Number: 107949GR
LOT Lot Number: NCUK220606 20524 001
Expiration Date: 2024-OCT-13

Reconstitute NHP Calibrator using 500 µL of Sample Diluent
Reconstitute NHP Custom Calibrator using 400 µL of NHP Calibrator

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	101852GR	NCUP240209	N/A	2025-AUG-09
Sample Diluent	101104HU	HQ6S230525	N/A	2025-MAY-25
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Calibrator 1	109255HU	NHPM220606	500 µL	2025-SEP-01
Calibrator 2	101855GR	NHPCUM221027	400 µL	2025-SEP-01
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	101856GR	NHPCUD240220	6 mL	2025-AUG-20
Streptavidin HRP	101173GR	QNSN231025	N/A	2025-APR-25
Substrate A	101159GR	QNSA231205	N/A	2024-OCT-13
Substrate B+	101120GR	QNSB+231205	N/A	2024-OCT-13
Wash Buffer	101158GR	QNSW230801	N/A	2025-AUG-01

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	CXCL-1	1760	pg/mL	1760	2.41	1.21
2	GMCSF	5970	pg/mL	5970	8.19	4.09
3	IFNγ	820	pg/mL	820	1.12	0.56
4	IL-1β	3660	pg/mL	3660	5.02	2.51
5	IL-2	560	pg/mL	560	0.77	0.38
6	IL-4	52	pg/mL	52	0.071	0.036
7	IL-5	300	pg/mL	300	0.41	0.21
8	IL-6	910	pg/mL	910	1.25	0.62
9	IL-8	1920	pg/mL	1920	2.63	1.32
10	IL-10	2460	pg/mL	2460	3.37	1.69
11	IL-12p70	2630	pg/mL	2630	3.61	1.80
12	IL-13	280	pg/mL	280	0.38	0.19
13	IL-15	3100	pg/mL	3100	4.25	2.13
14	IL-17	6390	pg/mL	6390	8.77	4.38
15	IL-18	186	pg/mL	186	0.26	0.13
16	TNFα	190	pg/mL	190	0.26	0.13
17	sIL-1RA	6490	pg/mL	6490	8.90	4.45
18	MCP-1	1880	pg/mL	1880	7.74	3.87
19	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: [Signature]

Date: 01 MAR 24

QA Release: [Signature]

Date: 01 MAR 24

for Research Use Only