


101896GR - 7.0

QUANSYS
BIOSCIENCESQuansys Biosciences
365 N 600 W Logan, UT 84321
(435) 752-0531**CERTIFICATE OF ANALYSIS**Q-Plex™ Porcine Cytokine HS (4-Plex)
Software Product Code: PCSM170306CY**REF** Part Number: 119149PC
LOT Lot Number: PCSK170306 20555 001
Expiration Date: 2025-MAR-13

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	119152PC	PCSP230609	N/A	2025-JUN-09
Sample Diluent	119157PC	PCSS230613	N/A	2025-JUN-13
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	N/A	N/A	N/A	N/A
Calibrator	119155PC	PCSM170306	500 µL	2025-JUN-21
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	119156PC	PCSD230615	N/A	2025-JUN-15
Streptavidin HRP	101173GR	QNSN240103	N/A	2025-JUL-03
Substrate A	101159GR	QNSA240111	N/A	2025-MAR-13
Substrate B+	101120GR	QNSB+240111	N/A	2025-MAR-13
Wash Buffer	101158GR	QNSW230314	N/A	2025-MAR-14

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	IL-1β	1600	pg/mL	N/A	N/A	2.10
2	IL-6	2000	pg/mL	N/A	N/A	2.50
3	IL-8	2300	pg/mL	N/A	N/A	2.90
4	TNFα	1350	pg/mL	N/A	N/A	1.90
5	REF	N/A	N/A	N/A	N/A	N/A

**Q-Plex™** ARRAYS

QC Release:

Date: 12MAR24

QA Release:

Date: 12MAR24

for Research Use Only