

101896GR - 7.0

QUANSYS 
 BIOSCIENCES
CERTIFICATE OF ANALYSIS


Quansys Biosciences
 365 N 600 W Logan, UT 84321
 (435) 752-0531

Q-Plex™ Human Custom 3-Plex
 Software Product Code: HA1M220712-L3 


REF Part Number: 107749GR-03
LOT Lot Number: HCUK220712 20991 001
 Expiration Date: 2025-FEB-02

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	101852GR	HCUP240625	N/A	2026-JUN-25
Sample Diluent	101104HU	HQ6S240125	N/A	2026-JAN-25
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	N/A	N/A	N/A	N/A
Calibrator	102655HU	HA1M220712	600 µL	2026-JUN-27
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	110956HU 01	HCSDQ6230612	N/A	2025-JUN-12
Streptavidin HRP	101173GR	QNSN240103	N/A	2025-JUL-03
Substrate A	101159GR	QNSA231027	N/A	2025-AUG-23
Substrate B+	101120GR	QNSB+231027	N/A	2025-AUG-23
Wash Buffer	101158GR	QNSW230202	N/A	2025-FEB-02


Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	IL-1β	11800	pg/mL	11800	16.19	8.09
2	IL-6	2850	pg/mL	2850	3.91	1.95
3	TNFα	6020	pg/mL	6020	8.26	4.13
4	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: 

Date: 27 JUN 24

QA Release: 

Date: 27 JUN 24

for Research Use Only