

101896GR - 7.0

QUANSYS
 BIOSCIENCES
CERTIFICATE OF ANALYSIS


Quansys Biosciences
 365 N 600 W Logan, UT 84321
 (435) 752-0531

Q-Plex™ Human Pneumococcal IgG (9-Plex)
 Software Product Code: HPNM231211-9P

REF Part Number: 481649HU
LOT Lot Number: HPNK231211 21072 001
 Expiration Date: 2024-NOV-28

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	481652HU	HPNP240618	N/A	2025-JUN-18
Sample Diluent	481157HU	HPNS240401	N/A	2025-APR-01
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	N/A	N/A	N/A	N/A
Calibrator	481155HU	HPNM231211	500 µL	2025-JUL-23
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	481156HU	HPND231128	N/A	2024-NOV-28
Streptavidin HRP	N/A	N/A	N/A	N/A
Substrate A	101159GR	QNSA240405	N/A	2026-MAR-19
Substrate B+	101120GR	QNSB+240405	N/A	2026-MAR-19
Wash Buffer	101158GR	QNSW240412	N/A	2026-APR-12

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	2	2.25	µg/mL	0.56	0.0022	0.0011
2	10A	1.19	µg/mL	0.070	0.0012	0.00058
3	11A	0.47	µg/mL	0.12	0.00046	0.00023
4	15B	1.55	µg/mL	0.39	0.0015	0.00076
5	17F	0.73	µg/mL	0.18	0.0028	0.0014
6	19A	1.28	µg/mL	0.32	0.0050	0.0025
7	20	0.92	µg/mL	0.23	0.00090	0.00045
8	22F	0.87	µg/mL	0.22	0.0034	0.0017
9	33F	0.94	µg/mL	0.060	0.00092	0.00046
10	CWPS	N/A	N/A	N/A	N/A	N/A
11	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release:

[Signature]

Date:

25 JUL 24

QA Release:

[Signature]

Date:

25 JUL 24

for Research Use Only