

101896GR - 7.0



Quansys Biosciences
365 N 600 W Logan, UT 84321
(435) 752-0531

CERTIFICATE OF ANALYSIS

Q-Plex™ Mouse Cytokine Panel 1 HS (6-Plex)

Software Product Code: MHSM230927-P16P


REF Part Number: 111349MS

LOT Lot Number: MCP1K230927 21060 001

Expiration Date: 2024-OCT-13

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	111352MS	MCP1P240212	N/A	2025-FEB-12
Sample Diluent	112457MS	MHSS240213	N/A	2025-JUN-06
Calibrator Diluent	N/A	N/A	N/A	N/A
Assay Diluent	N/A	N/A	N/A	N/A
Competitor	N/A	N/A	N/A	N/A
Calibrator	112455MS	MHSM230927	550 µL	2025-JUL-19
Calibrator Additive	N/A	N/A	N/A	N/A
Controls	N/A	N/A	N/A	N/A
Detection	112456MS	MHSD231013	6 mL	2024-OCT-13
Streptavidin HRP	101173GR	QNSN240508	N/A	2025-SEP-06
Substrate A	101159GR	QNSA240405	N/A	2026-MAR-19
Substrate B+	101120GR	QNSB+240405	N/A	2026-MAR-19
Wash Buffer	101158GR	QNSW230202	N/A	2025-FEB-02

Spot #	Analyte	Cal Conc	Units	ULOQ	LLOQ	LLD
1	IFN γ	800	pg/mL	800	1.10	0.55
2	IL-1 α	620	pg/mL	620	0.85	0.43
3	IL-1 β	2390	pg/mL	2390	3.28	1.64
4	IL-6	1110	pg/mL	1110	1.52	0.76
5	IL-10	520	pg/mL	520	0.71	0.36
6	TNF α	390	pg/mL	390	0.53	0.27
7	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release: Date: 19 JUL 24
 QA Release: Date: 19 JUL 24

for Research Use Only