

101896GR - 7.0



Quansys Biosciences
365 N 600 W Logan, UT 84321
(435) 752-0531


CERTIFICATE OF ANALYSIS

Q-Plex™ Human Pneumococcal IgG (9-Plex) Aliquoted Bulk 30 mL Diluent
Software Product Code: HPNM240514-A9P

REF Part Number: 481249HU-25
LOT Lot Number: HPNK240514 21168 001
Expiration Date: 2025-JUN-11

| Component | Part Number | Lot Number | Recon. Volume | Exp. Date |
|---------------------|-------------|-------------|---------------|-------------|
| Plate | 481652HU | HPNP240618 | N/A | 2025-JUN-18 |
| Sample Diluent | 481157HU-30 | HPNS240611 | N/A | 2025-JUN-11 |
| Calibrator Diluent | N/A | N/A | N/A | N/A |
| Assay Diluent | N/A | N/A | N/A | N/A |
| Competitor | N/A | N/A | N/A | N/A |
| Calibrator | 481155HU | HPNM240514 | 500 µL | 2025-AUG-20 |
| Calibrator Additive | N/A | N/A | N/A | N/A |
| Controls | N/A | N/A | N/A | N/A |
| Detection | 481156HU | HPND240621 | N/A | 2025-JUN-21 |
| Streptavidin HRP | N/A | N/A | N/A | N/A |
| Substrate A | 101159GR | QNSA240531 | N/A | 2026-MAR-19 |
| Substrate B+ | 101120GR | QNSB+240531 | N/A | 2026-MAR-19 |
| Wash Buffer | 101158GR | QNSW240809 | N/A | 2026-AUG-09 |

| Spot # | Analyte | Cal Conc | Units | ULOQ | LLOQ | LLD |
|--------|---------|----------|-------|-------|---------|---------|
| 1 | 2 | 2.19 | µg/mL | 0.43 | 0.0031 | 0.0016 |
| 2 | 10A | 1.16 | µg/mL | 0.057 | 0.0017 | 0.00085 |
| 3 | 11A | 0.44 | µg/mL | 0.091 | 0.00060 | 0.00030 |
| 4 | 15B | 1.51 | µg/mL | 0.29 | 0.0021 | 0.0011 |
| 5 | 17F | 0.74 | µg/mL | 0.14 | 0.0043 | 0.0022 |
| 6 | 19A | 1.36 | µg/mL | 0.25 | 0.0073 | 0.0037 |
| 7 | 20 | 0.91 | µg/mL | 0.18 | 0.0014 | 0.00070 |
| 8 | 22F | 0.85 | µg/mL | 0.17 | 0.0048 | 0.0024 |
| 9 | 33F | 0.96 | µg/mL | 0.047 | 0.0014 | 0.00070 |
| 10 | CWPS | N/A | N/A | N/A | N/A | N/A |
| 11 | REF | N/A | N/A | N/A | N/A | N/A |



Q-Plex™ ARRAYS

QC Release: [Signature]

Date: 20 AUG 24

QA Release: [Signature]

Date: 20 AUG 24

for Research Use Only