



101896GR - 8.0

QUANSYS
 BIOSCIENCES
CERTIFICATE OF ANALYSIS


Quansys Biosciences
 365 N 600 W
 Logan, UT 84321
 (435) 752-0531

Q-Plex™ Human T Cell Cytokine HS (9-Plex)

REF	Part Number: 111933HU	Software Product Code: HCHM200803-T9 
LOT	Lot Number: HTCK200803 21324 001	
	Expiration Date: 2025-SEP-06	

Reconstitute Calibrator using 750 µL of Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	111952HU	HTCP240619	N/A	2026-JUN-19
Sample Diluent	101104HU	HQ6S240125	N/A	2026-JAN-25
Calibrator	112455HU	HCHM200803	750 µL	2026-SEP-19
Detection	101856GR	HCHDQ6240614	6 mL	2026-JUN-14
Streptavidin HRP	101173GR	QNSN240508	N/A	2025-SEP-06
Substrate A	101159GR	QNSA240405	N/A	2026-MAR-19
Substrate B+	101120GR	QNSB+240405	N/A	2026-MAR-19
Wash Buffer	101158GR	QNSW240702	N/A	2026-JUL-02



Q-Plex™ ARRAYS

Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	IFN γ	93	pg/mL	0.064	0.13	93
2	IL-2	47	pg/mL	0.032	0.064	47
3	IL-4	18	pg/mL	0.012	0.025	18
4	IL-5	260	pg/mL	0.18	0.36	260
5	IL-10	2320	pg/mL	1.59	3.18	2320
6	IL-13	290	pg/mL	0.20	0.40	290
7	IL-17	4500	pg/mL	3.09	6.17	4500
8	IL-21	260	pg/mL	0.18	0.36	260
9	TNF α	640	pg/mL	0.44	0.88	640
10	REF	N/A	N/A	N/A	N/A	N/A

QC Release:

J. Wa

Date:

19 SEP 24

QA Release:

Do

Date:

19 SEP 24

for Research Use Only