


101896GR - 8.0

QUANSYS
BIOSCIENCES
CERTIFICATE OF ANALYSISQuansys Biosciences
365 N 600 W
Logan, UT 84321
(435) 752-0531**Q-Plex™ Human Cytokine Inflammation Panel 1 (9-Plex)**

REF	Part Number: 110433HU	Software Product Code: HA1M220712-IF 
LOT	Lot Number: HClK220712 21274 001	
	Expiration Date: 2025-JUN-06	

Reconstitute Calibrator using 600 µL of Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	110452HU	HCIP230606	N/A	2025-JUN-06
Sample Diluent	101104HU	HQ6S240125	N/A	2026-JAN-25
Calibrator	102655HU	HA1M220712	600 µL	2026-SEP-13
Detection	110956HU 01	HCSDQ6230612	N/A	2025-JUN-12
Streptavidin HRP	101173GR	QNSN240508	N/A	2025-SEP-06
Substrate A	101159GR	QNSA240531	N/A	2026-MAR-19
Substrate B+	101120GR	QNSB+240531	N/A	2026-MAR-19
Wash Buffer	101158GR	QNSW240702	N/A	2026-JUL-02



Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	IL-1α	3940	pg/mL	3.55	7.10	2800
2	IL-1β	11800	pg/mL	10.90	21.80	8558
3	IL-2	5150	pg/mL	4.65	9.30	3649
4	IL-4	1670	pg/mL	1.50	3.00	1173
5	IL-6	2850	pg/mL	2.65	5.30	2100
6	IL-8	2080	pg/mL	1.80	3.60	1400
7	IL-10	8400	pg/mL	7.15	14.30	5627
8	IFNγ	9775	pg/mL	7.15	14.30	5600
9	TNFα	6020	pg/mL	5.75	11.50	4522
10	REF	N/A	N/A	N/A	N/A	N/A

Q-Plex™ ARRAYS

QC Release:

QA Release:

Date:

Date:

for Research Use Only

16 SEP 24

16 Sep 24