



101896GR - 8.0

QUANSYS BIOSCIENCES CERTIFICATE OF ANALYSIS

Quansys Biosciences
365 N 600 W
Logan, UT 84321
(435) 752 - 0531


Q-Plex™ Human Custom Singleplex

REF	Part Number: 107749GR-01	Software Product Code: HA1M220719-IL6 
LOT	Lot Number: HCUK220719 21519 001	
	Expiration Date: 2026-JAN-01	

Reconstitute Calibrator using 650 µL of Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	101852GR	HCUP241024	N/A	2026-OCT-24
Sample Diluent	101104HU	HQ6S240805	N/A	2026-AUG-05
Calibrator	102655HU	HA1M220719	650 µL	2026-NOV-01
Detection	101856GR	HCUD240628	6 mL	2026-JUN-28
Streptavidin HRP	101173GR	QNSN240701	N/A	2026-JAN-01
Substrate A	101159GR	QNSA240801	N/A	2026-MAR-19
Substrate B+	101120GR	QNSB+240801	N/A	2026-MAR-19
Wash Buffer	101158GR	QNSW240809	N/A	2026-AUG-09

Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	IL-6	2710	pg/mL	1.86	3.71	2710
2	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release:

P. Wa

Date:

01 NOV 24

QA Release:

Stacy Hansen
for Research Use Only

Date:

01 NOV 24