




Q-Plex™ Human T Cell Cytokine HS (9-Plex)

REF	Part Number: 111933HU	Software Product Code: HCHM200803-T9 
LOT	Lot Number: HTCK200803 21619 001	
	Expiration Date: 2025-MAR-26	

Reconstitute Calibrator using 750 µL of Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	111952HU	HTCP240619	N/A	2026-JUN-19
Sample Diluent	101104HU	HQ6S240805	N/A	2026-AUG-05
Calibrator	112455HU	HCHM200803	750 µL	2026-NOV-25
Detection	112456HU 01	HCHDQ6230326	6 mL	2025-MAR-26
Streptavidin HRP	101173GR	QNSN241105	N/A	2026-MAY-05
Substrate A	101159GR	QNSA240801	N/A	2026-MAR-19
Substrate B+	101120GR	QNSB+240801	N/A	2026-MAR-19
Wash Buffer	101158GR	QNSW240809	N/A	2026-AUG-09



Q-Plex™ ARRAYS

Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	IFN γ	93	pg/mL	0.064	0.13	93
2	IL-2	47	pg/mL	0.032	0.064	47
3	IL-4	18	pg/mL	0.012	0.025	18
4	IL-5	260	pg/mL	0.18	0.36	260
5	IL-10	2320	pg/mL	1.59	3.18	2320
6	IL-13	290	pg/mL	0.20	0.40	290
7	IL-17	4500	pg/mL	3.09	6.17	4500
8	IL-21	260	pg/mL	0.18	1.00	260
9	TNF α	640	pg/mL	0.44	0.88	640
10	REF	N/A	N/A	N/A	N/A	N/A

QC Release: *J. Wa*

Date: 25 NOV 24

QA Release: *K. ax*

Date: 25 NOV 24

for Research Use Only