



101896GR - 8.0

# QUANSYS BIOSCIENCES CERTIFICATE OF ANALYSIS


Quansys Biosciences  
365 N 600 W  
Logan, UT 84321  
(435) 752-0531

## Q-Plex™ Mouse Cytokine Panel 3 HS (4-Plex)

<b>REF</b>	Part Number: 111249MS	Software Product Code: <b>MHSM230927-P34P</b> 
<b>LOT</b>	Lot Number: MCP3K230927 21741 001	
	Expiration Date: 2025-SEP-19	

### Reconstitute Calibrator using 550 µL of Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	111252MS	MCP3P250107	N/A	2026-JAN-07
Sample Diluent	112457MS	MHSS240730	N/A	2026-JUL-30
Calibrator	112455MS	MHSM230927	550 µL	2026-JAN-13
Detection	112456MS	MHSD240919	6 mL	2025-SEP-19
Streptavidin HRP	101173GR	QNSN241105	N/A	2026-MAY-05
Substrate A	101159GR	QNSA240801	N/A	2026-MAR-19
Substrate B+	101120GR	QNSB+240801	N/A	2026-MAR-19
Wash Buffer	101158GR	QNSW240809	N/A	2026-AUG-09



Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	IL-2	180	pg/mL	0.12	0.25	180
2	IL-4	160	pg/mL	0.11	0.22	160
3	IL-6	1110	pg/mL	0.76	1.52	1110
4	TNFα	390	pg/mL	0.27	0.53	390
5	REF	N/A	N/A	N/A	N/A	N/A

Q-Plex™ ARRAYS

QC Release:

Date:

13 JAN 25

QA Release:

Date:

13 JAN 25

for Research Use Only