

101896GR - 8.1

QUANSYS
BIOSCIENCES
CERTIFICATE OF ANALYSIS


Quansys Biosciences
365 N 600 W
Logan, UT 84321
(435) 752 - 0531

Q-Plex™ Human Custom 5-Plex

REF	Part Number: 107749GR-05	Software Product Code: HA1M220719-L5 
LOT	Lot Number: HCUK220719 22096 001	
	Expiration Date: 2026-MAR-19	

Reconstitute Calibrator 1 using 650 µL of Sample Diluent
Reconstitute Calibrator 2 using 400 µL of Calibrator 1

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	101852GR	HCUP250225	N/A	2027-FEB-25
Sample Diluent	101104HU	HQ6S240805	N/A	2026-AUG-05
Calibrator 1	102655HU	HA1M220719	650 µL	2027-MAR-27
Calibrator 2	102855HU	HA2M200303	400 µL	2027-MAR-27
Detection	101856GR	HCUD250321	6 mL	2027-MAR-21
Streptavidin HRP	101173GR	QNSN241105	N/A	2026-MAY-05
Substrate A	101159GR	QNSA241108	N/A	2026-MAR-19
Substrate B+	101120GR	QNSB+241108	N/A	2026-MAR-19
Wash Buffer	101158GR	QNSW250106	N/A	2027-JAN-06



Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	IFN α	340	pg/mL	0.23	0.47	340
2	IFN β	21700	pg/mL	14.88	29.77	21700
3	IFN γ	8520	pg/mL	5.84	11.69	8520
4	IFN λ	12480	pg/mL	8.56	17.12	12480
5	IFN ω	3990	pg/mL	2.74	5.47	3990
6	REF	N/A	N/A	N/A	N/A	N/A

Q-Plex™ ARRAYS

QC Release:

Date: 27MAR25

QA Release:

Date: 27MAR25

for Research Use Only